



Getting Started
Grid Pad Go 8 and 10

Looking after your Grid Pad

When not in use, store your Grid Pad safely with no hard objects rested or pressed against the screen.

Cleaning

Keeping your device clean is important, particularly if it is being shared among multiple users or used in a school or hospital.

You can wipe the screen with a damp cloth or antibacterial wipe. Tough grime can be dislodged with an old toothbrush or similar.

To clean the ports, you can use a toothbrush in combination with a vacuum cleaner to remove any dirt. Do not make the ports wet. Be careful using a vacuum cleaner around any vents on your device, this can damage the internals.


Do not pour water on or submerge your communication aid!

Protection

The Grid Pad Go is rated to IP-65 for protection from splashes and drops, with sealed ports and tough outer case. The device is built to withstand a drop of 1.2m.


Turning on your Grid Pad

To turn on your Grid Pad press the power button on the top of the device:



Charging your Grid Pad

To charge your Grid Pad, plug the charging cable into a power source and insert into the device as shown below.



Pairing the Bluetooth speaker

If you need to re-pair your Bluetooth speaker with your Grid Pad please follow the steps below:

1. Follow the instructions provided with your speaker to enter pairing mode.
2. On your Grid Pad go to:
Windows Settings > PC and Devices > Bluetooth > Bluetooth On. It will now automatically search.
3. Select your speaker from the list to pair.

The speaker is charged by connecting the USB port to Grid Pad PSU.

Battery guidance

Getting the most out of your battery

The Grid Pad Go battery will last an average user 7 hours. Every user will use the device in a different way. Intensive use including using Bluetooth, streaming video and listening to music for long periods of time will use the battery faster.

Charging guidance

Your devices battery will take between 2.5 to 3.5 hours to charge.

Specification

Go 10

Processor	Intel Atom E3940 Quad Core 1.6 GHz (1.8 GHz Turbo)
Memory	4GB
Hard drive	64GB
Graphics card	HD 500 Base 400 MHz, Burst 600 MHz
Dimensions (mm)	269 x 180 x 12.5
Battery	10 hours
Camera	2 MP front and 8 MP rear
Screen	10.1" 2560 x 1600
Connectivity	Wi-fi and Bluetooth
Ports	1 x USB-C

Specification

Go 8

Processor	Intel Atom E3940 Quad Core 1.6 GHz (1.8 GHz Turbo)
Memory	4GB
Hard drive	64GB
Graphics card	HD 500 Base 400 MHz, Burst 600 MHz
Dimensions (mm)	228 x 150 x 13
Battery life	7 hours
Camera	8 MP rear
Screen	8.4" - 2560 x 1600
Connectivity	Wi-fi and Bluetooth
Ports	1 x USB-C

Warranty

Grid Pad devices are covered under the standard warranty for 2 years.

Repairs

Please do not attempt to open your Grid Pad, as there are no user serviceable parts inside. If you want to discuss a repair please contact your local reseller.

For more information visit www.thinksmartbox.com

Statement of compliance

This apparatus complies with the requirements of directive 93/42/EEC concerning medical devices. This Product is subject to Directive 2002/96/EC of the European Parliament and the Council of the European Union on waste electrical and electronic equipment (WEEE), and should not be disposed of as unsorted municipal waste. Please utilise your local WEEE collection facilities in the disposal of this product.



Smartbox

US Office

Smartbox Assistive Technology Inc
2831 Leechburg Road
Lower Burrell
Westmoreland
PA
15068-2539

UK Office

Smartbox Assistive Technology Ltd
Ysobel House, Enigma Commercial Centre
Sandys road
Malvern, WR14 1JJ, UK

+44 (0) 1684 578868
info@thinksmartbox.com

www.thinksmartbox.com