

Manual

Created by

Smartbox

This guide provides the basics for getting started with your Touch Pad.



hub.thinksmartbox.com/TouchPad

### **Smart**box

Document: Touch Pad - Manual - UK - K

Language: English

Date of issue: Decmeber 2023

Copyright © 2023 Smartbox Assistive Technology



### Manual

1.	Starting Grid 3	4
2.	Power	6
3.	Device guide	8
4.	Charging	10
5.	Mounting	12
6.	Shoulder strap	14

# 1 Starting Grid 3

Grid 3 will automatically start when the device is switched on by default. However, you can start Grid 3 from icon on the desktop.



### For more resources and guides on how to get started with Grid visit the Smartbox Hub.


- Get started with Grid
- Access methods and setup
- Troubleshooting guides
- Videos and webinars
- Resources for grid sets
- ...and much more.



hub.thinksmartbox.com/grid3


# 2 Power

Before switching on your Touch Pad for the first time, it is best to plug in the power lead and charge the battery fully. Press the power button to turn the device on.



The device will boot into Windows, and Grid 3 will automatically start.

While the device boots, also press the Speaker power button to turn on the bluetooth speaker.



When not in use, it is good practice to shut down your Touch Pad completely. This will give any automatic updates a chance to be applied. We also recommend that your device is shut down before it is stored for transit


Open the Start menu and tap the power icon. From here you can select Sleep, Shut down and Restart. Selecting Shut down will power the device down completely.


You can also turn off your device from within Grid 3, by using the Turn off computer command. This command can be added to any grid set.

#### Force shut down

If you press and hold the power button for ten seconds, your Touch Pad will shut down. Please note that this is only recommended when the Touch Pad is non-responsive. Any unsaved work will be lost.








# 4 Charging

To charge your Touch Pad, plug the power lead into the charging port on the left side of the device.

To charge both the Touch Pad and the speaker at the same time use the provided charging lead with USB C.

A full charge will take up to 3 hours.



#### Using a desk, floor, or wheelchair mount


To mount your Touch Pad, you will need the optional dual mounting plate that will attach to either Daessy or ReHadapt.

#### Attaching the mounting plate

To attach the mounting plate, lay your Touch Pad with the screen facing down on a flat surface. To prevent damage to the screen, place the device on something soft like a towel, cloth or newspaper. You will need the supplied Allen key, mounting plate and the three screws.

Remove the stickers covering the three screw holes. Use the supplied Allen key to remove the three screws located on the back of the Touch pad, as shown below.





### Mounting your Touch Pad on a desk or floor stand

With the correct mounting plate attached, you can fit the Touch Pad to your wheelchair mount, floor stand or desk stand. As you attach the device to the stand, be sure to hold the Touch Pad securely until the head of the stand is firmly docked with the mounting plate.



For further instructions on positioning see the documentation provided with your stand.


### Shoulder strap

Your Touch Pad is supplied with a shoulder strap. These clip to the rings on the back of the device and can be easily attached and removed as you need.

First loop the clips to the case of the Touch Pad.

Then you can attach the longer strap.

The strap can be adjusted to a suitable length.







hub.thinksmartbox.com/TouchPad

**Smart**box

thinksmartbox.com

Touch Pad - Manual - UK - K