

LOOK Lab

Workbook

Created by

Smartbox

Introduction

Look Lab is a collection of activities to introduce and improve a range of eye gaze skills. More information about Look Lab can be found in the **Getting Started Guide**. This workbook is designed to be used alongside the software.

Worksheets

This workbook contains a worksheet to use alongside the activities. You can record progress, set objectives and assess engagement, independence and success in each session.

At the end of the workbook is a General record sheet, which is designed to act as a catch all sheet for any adhoc observations.

Easy and Standard modes

Most activities have the option to choose **Easy** or **Standard** mode.

Easy mode: games are usually slower, with larger target areas. Speed of play remains constant. We have also removed timers and lives.

Standard mode: some games are set to speed up and become harder as the play progresses. Games may have limited time or lives.



On each worksheet, we have listed the difference between the two modes and the effect it has on the activity.

Minimal and Detailed visual modes

You can choose between **Minimal** or **Detailed** mode via the **Options** button, in the Accessibility section.

Minimal mode: less distractions, with plain coloured backgrounds. The contrast is higher, making objects on screen easier to see.

Detailed mode: colourful backgrounds and images with the running score often displayed throughout the game.

Analysis

You can access a '**heat map**' via the Analysis button when an activity is paused or ended. This shows which areas of the screen have been looked at, and for how long.

These can be used to evidence access to different parts of the screen, and responses to screen activity or verbal prompts. They are a valuable tool which can be used alongside others (visual observations, recordings, scores etc) to track progress.

Contents

Skills list

4



Notes

Skills Checklists

Fun

6

Challenge

7

Puzzle

8

Mindfulness

9

Worksheets

Fun

11 - 21

Challenge

22 - 31

Puzzle

32 - 37

Mindfulness

38 - 43

Creative

44 - 48

Quick record sheet

49

Skills list

The following list explains the range of skills which Look Lab can help to develop when learning to use eye gaze:



1

Exploration

These activities encourage the user to look around the entire screen

2

Cause & Effect

These activities allow the user to see an immediate effect as they move their gaze around the screen

3

Tracking

These activities require the user to look at one or more objects on the screen and follow them with their eyes as they move.

4

Sustained Attention

These activities require the user to focus on certain areas of the screen with no distractions from other areas.

5

Dwelling

These activities require the user to fix their gaze on a target to make a selection.

6

Choice Making

These activities allow the user to choose between options according to their preference.

7

Reactions

These activities challenge the user to respond to sudden changes on screen.

8

Memory

These activities require the user to remember where objects are.

Notes



9 Horizontal Control

These activities require the user to move their gaze from left to right, encouraging the user to look at, or move, objects horizontally across the screen.

10 Vertical Control

These activities require the user to move their gaze up and down only. These encourage the user to look at, or move, objects vertically across the screen.

11 Precision

These activities require the user to display a high level of precise eye gaze skills to achieve success.

12 Problem Solving

These activities require logical thinking to overcome the challenge and make planned eye gaze movements to solve puzzles.

Skills Checklist



Fun: relaxing games that are easy to learn and get started with.

Activity	Exploration	Cause & Effect	Tracking	Sustained attention	Dwelling	Choice making	Reactions	Memory	Horizontal control	Vertical control	Precision	Problem solving	Compensation
Breakthrough			✓				✓			✓			
Starfish			✓	✓			✓						
Circle tennis			✓				✓						
Tunnel Dash				✓			✓				✓		
Hoops			✓				✓		✓				
Block Hunter				✓	✓							✓	
Paddle Ball			✓	✓									
Star Gaze	✓	✓	✓										
Penalty Spot			✓				✓		✓				
Sweet Slice!		✓	✓				✓						
Pocket Ball	✓										✓	✓	

Skills Checklist



Challenge: more complex and strategic games that involve time dependent reactions.

Activity	Exploration	Cause & Effect	Tracking	Sustained attention	Dwelling	Choice making	Reactions	Memory	Horizontal control	Vertical control	Precision	Problem solving	Compensation
Pipe Dodge				✓			✓			✓			
Bubble Jump							✓				✓	✓	
Crawly Cruncher	✓			✓			✓						
Shape Shifter							✓			✓		✓	
Space Shoot			✓		✓						✓		
Leap Frog					✓		✓					✓	
Zombie Splash!				✓	✓		✓						
Aliens!			✓				✓		✓				
Hedge Hopper					✓	✓						✓	
Splat	✓				✓		✓				✓	✓	

Skills Checklist



Puzzle: activities with a problem-solving element.

Activity	Exploration	Cause & Effect	Tracking	Sustained attention	Dwelling	Choice making	Reactions	Memory	Horizontal control	Vertical control	Precision	Problem solving	Compensation
Match n Smash					✓						✓	✓	
Memory Match					✓			✓					
Slide Puzzle					✓							✓	
Shuffle cups			✓	✓	✓								
Simon Says					✓			✓					
Mouse Maze					✓							✓	

Skills Checklist



Mindfulness: simple cause and effect activities that support relaxation and screen engagement.

Activity	Exploration	Cause & Effect	Tracking	Sustained attention	Dwelling	Choice making	Reactions	Memory	Horizontal control	Vertical control	Precision	Problem solving	Compensation
Fish Pond	✓	✓			✓								
Fireworks	✓	✓			✓								
Popcorn	✓	✓			✓								
Ball Pit	✓	✓			✓								
Pattern Pop	✓	✓									✓		
Chimes		✓							✓				

Skills Checklist



Creative: activities focused on music or art.

Activity	Exploration	Cause & Effect	Tracking	Sustained attention	Dwelling	Choice making	Reactions	Memory	Horizontal control	Vertical control	Precision	Problem solving	Compensation
Beatbox					✓	✓							
Guitar					✓	✓					✓		
Piano					✓	✓					✓		
Pixel Paint	✓										✓	✓	
Colour In	✓					✓					✓		



Star Gaze

Relax and enjoy a galaxy of dancing stars. Catch them all to reveal a new pattern!

Skills: Exploration, Cause and effect, Tracking

In Easy mode: The stars are bigger and move more slowly.



Skill	Objective	Not Observed	Developing	Achieved	Comments
Exploration	Moves gaze around the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Cause and effect	Looks at stars to cause a reaction	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Tracking	Follows the stars as they move around the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Hoops

Catch as many balls as you can but avoid the bombs!

Skills: Tracking, Reactions, Horizontal control

In Easy mode: The area in which the balls drop is narrower and there are no bombs.



Skill	Objective	Not Observed	Developing	Achieved	Comments
Tracking	Follows the balls as they appear on the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Moves the basket in time to catch the ball	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Horizontal control	Moves the basket side to side across the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Block Hunter

Find all the blocks which match the target and pop them!

Skills: Sustained attention, Dwelling, Problem Solving

In Easy mode: The blocks are bigger and there are fewer of them.



Skill	Objective	Not Observed	Developing	Achieved	Comments
Sustained attention	Focuses on a single block and is not distracted by other blocks moving	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Dwelling	Looks at a block for long enough to select it	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Problem solving	Selects blocks from the correct category only	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Breakthrough

Use your bat to bounce the ball and smash the blocks to clear the level. Don't miss the ball!

Skills: Tracking, Reactions, Vertical Control

In Easy mode: The ball moves more slowly, and the bat is bigger.



Fun

Skill	Objective	Not Observed	Developing	Achieved	Comments
Tracking	Follows the ball as it moves	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Moves the bat in time to receive the ball	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Vertical control	Moves the bat up and down the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Circle Tennis

Keep the ball inside the circle for as long as you can!

Skills: Tracking, Reactions

In Easy mode: The ball moves more slowly and the bat is bigger.



Skill	Objective	Not Observed	Developing	Achieved	Comments
Tracking	Follows the ball as it moves	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Moves the bat around as the ball changes direction	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Starfish

Stop the jellyfish from ruining the starfish's holiday. Collect stars to deliver a burst of positive energy!

Skills: Tracking, Sustained attention, Reactions

In Easy mode: The jellyfish move more slowly, and appear less often.



Fun

Skill	Objective	Not Observed	Developing	Achieved	Comments
Tracking	Follows the jellyfish around the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Sustained attention	Focuses on the jellyfish that is closest to the starfish	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Moves shield in time to repel the jellyfish	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Paddle Pop

Keep bouncing the ball on your paddle - don't let it drop

Skills: Tracking, Sustained attention

In Easy mode: The paddle is larger.



Fun

Skill	Objective	Not Observed	Developing	Achieved	Comments
Tracking	Follows the ball around the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	Focuses on the ball throughout the game	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Tunnel Dash

Guide the ball through the gates and collect as many stars as you can.

Skills: Sustained attention, Reactions, Precision

In Easy mode: The gates move more slowly and the gap in each one is bigger.



Skill	Objective	Not Observed	Developing	Achieved	Comments
Sustained attention	Focuses on the gate until they move through it	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Moves the ball in time to pick up stars and go through gates	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Precision	Moves the ball accurately through gates	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Sweet Slice

Slice through the cakes as they fly through the air. Don't let them fall.

Skills: Cause and effect, Tracking, Reactions

In Easy mode: The cakes move more slowly. If you miss, you can just carry on.



Fun

Skill	Objective	Not Observed	Developing	Achieved	Comments
Cause and effect	Looks at a cake to slice it	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Tracking	Follows the cakes as they move around	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Looks at a new cake as it appears	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Penalty Spot

Keep the ball in the air for as long as possible!

Skills: Tracking, Reactions, Horizontal control

In Easy mode: The ball moves more slowly through the air and does not speed up.



Fun

Skill	Objective	Not Observed	Developing	Achieved	Comments
Tracking	Follows the ball as it moves around the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Moves the platform in time to catch the ball	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Horizontal control	Moves the platform side to side across the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Pocket Ball

Drive the coloured balls into the matching holes.

Skills: Exploration, Precision, Problem solving



In Easy mode: There is no time penalty if the ball goes in the wrong pocket.

Skill	Objective	Not Observed	Developing	Achieved	Comments
Exploration	Looks all around the screen to direct the digger	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Precision	Controls the digger to push the ball accurately	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Problem solving	Uses the digger to push the ball into the right hole	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Crawly Cruncher

Guide the caterpillar to crunch stars, but avoid the spikey bugs!

Skills: Exploration, Sustained attention, Reactions

In Easy mode: No nasty creatures to avoid and the caterpillar's speed stays constant.



Challenge

Skill	Objective	Not Observed	Developing	Achieved	Comments
Exploration	Looks at different areas of the screen to move around	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Sustained attention	Focuses on stars to collect them	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Moves away from the spikey bugs	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name


Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)

16



Splat

Worms are ruining your lawn! Bop them back down into their worm holes as fast as you can.

Skills: Exploration, Dwelling, Reactions



Challenge

In Easy mode: There is no time limit.

Skill	Objective	Not Observed	Developing	Achieved	Comments
Exploration	Looks at different areas of the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Dwelling	Looks at each worm long enough to bop it	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Responds to new worms popping up	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Comments (e.g. Device position, user position, calibration)

Facilitator

Time



Pipe Dodge

Guide the bird between the pipes and collect as many coins as you can!

Skills: Sustained attention, Reactions, Vertical control

In Easy mode: Coins and the gaps in between pipes are bigger. The bird moves more slowly.



Challenge

Skill	Objective	Not Observed	Developing	Achieved	Comments
Sustained attention	Focuses on the bird rather than the pipes	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Responds to coins and pipes as they appear	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Vertical control	Moves the bird up and down the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Aliens

Aliens are invading! Protect the Earth by zapping them all.

Skills: Tracking, Reactions, Horizontal control

In Easy mode: Aliens move more slowly and there is no lava to avoid.



Challenge

Skill	Objective	Not Observed	Developing	Achieved	Comments
Tracking	Follows the aliens as they move around the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Responds to changes in direction and fireballs	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Horizontal control	Moves the spaceship across the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Space Shoot

Shoot the asteroids before they hit your spaceship, and shoot the stars to collect points

Skills: Tracking, Dwelling, Sustained attention

In Easy mode: The asteroids move more slowly.



Challenge

Skill	Objective	Not Observed	Developing	Achieved	Comments
Tracking	Follows the asteroids as they move	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Dwelling	Looks at an asteroid long enough to shoot it	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Sustained attention	Keeps focus on closest asteroid as others appear	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Zombie Splash

Take aim and fire your water pistol to soak the zombies until they drop. Don't let them get too close!

Skills: Sustained attention, Dwelling, Reactions

In Easy mode: Zombies move more slowly and there are fewer of them.



Challenge

Skill	Objective	Not Observed	Developing	Achieved	Comments
Tracking	Follows the asteroids as they move	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Dwelling	Looks at an asteroid long enough to shoot it	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Sustained attention	Keeps focus on closest asteroid as others appear	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Shape Shifter

Fly through the gate which matches your spaceship's colour and shape, but look out for gates which change your spaceship!

Skills: Reactions, Vertical control, Problem solving



Challenge

In Easy mode: There are fewer gates to choose from and the spaceship moves more slowly.

Skill	Objective	Not Observed	Developing	Achieved	Comments
Reactions	Responds to new gates by moving the spaceship	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Vertical control	Moves the spaceship up and down the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Problem solving	Guides the spaceship through the matching gate	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Leap Frog

Jump around to dodge the monster and collect stars. Use the dwell marker to plan your moves!

Skills: Dwelling, Reactions, Problem solving

In Easy mode: The monster is slower and the focus points are bigger.



Challenge

Skill	Objective	Not Observed	Developing	Achieved	Comments
Dwelling	Looks at focus point long enough to move the frog	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Reactions	Moves towards stars as they appear and away from the monster	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Problem solving	Targets stars whilst avoiding monster	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Hedge Hopper

Foxes are causing havoc in the neighbourhood. Help the hedgehogs to defend it!

Skills: Dwelling, Choice making, Problem solving

In Easy mode: There are unlimited hedgehogs in each round.



Challenge

Skill	Objective	Not Observed	Developing	Achieved	Comments
Dwelling	Looks at target long enough to fire	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Choice making	Chooses where to place target each time	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Problem solving	Plans where to aim available hedgehogs	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Bubble Jump

Bounce off the platforms and help the bubble climb as high as it can!
Earn extra points by collecting coins

Skills: Reactions, Precision, Problem solving



Challenge

In Easy mode: Platforms are larger and there are fewer of them on screen at any time.

Skill	Objective	Not Observed	Developing	Achieved	Comments
Reactions	Looks at platforms in time to direct the bubble	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Precision	Accurately directs the bubble to a platform	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Problem solving	Chooses an appropriate platform to jump to	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Match n Smash

Match three or more tiles to clear them before they fill the screen. Clear them faster by matching more at once!

Skills: Reactions, Precision, Problem solving

In Easy mode: The tiles fill the screen more slowly and there are fewer shapes / colours.



Puzzle

Skill	Objective	Not Observed	Developing	Achieved	Comments
Dwelling	Looks at destination long enough to fire block	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Precision	Looks accurately at destination	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Problem solving	Places blocks to build up to three in a row	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Memory Match

Find and match the pairs of cards!

Skills: Dwelling, Memory

In Easy mode: There will always only be 4 cards.



Puzzle

Skill	Objective	Not Observed	Developing	Achieved	Comments
Dwelling	Focuses on a card long enough to turn it over	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Memory	Looks accurately at destination	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Slide Puzzle

Slide the puzzle pieces to reveal the jumbled up picture

Skills: Dwelling, Problem solving



Puzzle

No Easy / Standard mode for this activity.

Skill	Objective	Not Observed	Developing	Achieved	Comments
Dwelling	Looks at a piece long enough to make it slide	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Problem solving	Moves pieces to unscramble the image	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Shuffle Cups

Watch the cups carefully as they move to keep track of the ball.

Skills: Tracking, Dwelling, Sustained attention

In Easy mode: There is only one shuffle and you will get a hint if you choose the wrong cup.



Puzzle

Skill	Objective	Not Observed	Developing	Achieved	Comments
Tracking	Follows the cup with the ball under it	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Sustained attention	Keeps focus on target cup throughout the shuffle	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Dwelling	Focuses on selected cup long enough to lift it up	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Simon Says

Watch the pattern carefully, then repeat it. The pattern gets longer each time!

Skills: Dwelling, Memory



Puzzle

In Easy mode: If you get the sequence wrong, you can try again as many times as you need.

Skill	Objective	Not Observed	Developing	Achieved	Comments
Dwelling	Looks at selected button long enough to activate it	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Memory	Remembers the sequence correctly	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Comments (e.g. Device position, user position, calibration)



Mouse Maze

Grab the cheese to get the keys!
Plan your moves carefully to show the mouse where to go.

Skills: Dwelling, Problem Solving

No Easy / Standard mode for this activity.



Puzzle

Skill	Objective	Not Observed	Developing	Achieved	Comments
Dwelling	Focuses on arrow long enough to select	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Problem Solving	Plans movement of the mouse to get the cheese	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Fish Pond

Relax and enjoy a moment of tranquility watching the fish in the pond. You can even throw in some food to attract them!

Skills: Exploration, Cause and Effect, Dwelling



Mindfulness

No Easy / Standard mode for this activity.

Skill	Objective	Not Observed	Developing	Achieved	Comments
Exploration	Looks around the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Cause and Effect	Looks at a fish to make it move	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Dwelling	Fixes gaze in one place to drop food	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Fireworks

Time to celebrate! Guide the fireworks over the city and put on the greatest display!

Skills: Exploration, Cause and Effect, Dwelling



Mindfulness

No Easy / Standard mode for this activity.

Skill	Objective	Not Observed	Developing	Achieved	Comments
Exploration	Looks around the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Cause and Effect	Sets off multiple fireworks	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Dwelling	Fixes gaze in one place to set off firework	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Popcorn

Push the tumbling popcorn around the screen and enjoy the show!

Skills: Exploration, Cause and Effect



Mindfulness

No Easy / Standard mode for this activity.

Skill	Objective	Not Observed	Developing	Achieved	Comments
Exploration	Looks around the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Cause and Effect	Looks at a popcorn to make it move	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Ball Pit

Dive in and bounce the colourful balls around the ball pit!

Skills: Exploration, Cause and Effect

No Easy / Standard mode for this activity.



Mindfulness

Skill	Objective	Not Observed	Developing	Achieved	Comments
Exploration	Looks around the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Cause and Effect	Looks at balls to make them move	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Pattern Pop

Take a moment to enjoy popping every single bubble on the screen. Pop, pop, pop!

Skills: Exploration, Cause and Effect, Precision

In Easy mode: The bubbles are larger.



Mindfulness

Skill	Objective	Not Observed	Developing	Achieved	Comments
Exploration	Looks around the screen	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Cause and Effect	Looks at bubbles to make them pop	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Precision	Gazes accurately to pop all of the bubbles	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Chimes

Enjoy soft wooden melodies as the chimes sway gently in the wind.

Skills: Cause and Effect, Horizontal control



Mindfulness

No Easy / Standard mode for this activity.

Skill	Objective	Not Observed	Developing	Achieved	Comments
Cause and Effect	Looks at the chimes to play a sound	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Horizontal control	Looks side to side on the screen to move chimes	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Beat Box

Build your own beat and enjoy the groove!

Skills: Dwelling, Choice making

In Easy mode: The cells are larger and the tempo is slower.



Creative

Skill	Objective	Not Observed	Developing	Achieved	Comments
Dwelling	Looks at chosen square long enough to select	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Choice making	Chooses where to place the beats	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Guitar

Choose your chords and play a song! Strum up and down, or let the guitar play by itself.

Skills: Dwelling, Choice Making, Precision

In Easy mode: The interface is simplified by having a reduced number of options.




Creative

Skill	Objective	Not Observed	Developing	Achieved	Comments
Dwelling	Looks at letters long enough to play a chord	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Choice making	Chooses which chord to play	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Precision	Plays individual strings	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Comments (e.g. Device position, user position, calibration)



Piano

Play the piano, pick out notes, chords or arpeggios. Choose from five different instrument sounds.

Skills: Dwelling, Choice Making, Precision

In Easy mode: The interface is simplified by having a reduced number of options.



Skill	Objective	Not Observed	Developing	Achieved	Comments
Dwelling	Looks at letters long enough to play a note	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Choice making	Chooses different modes and instruments	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Precision	Looks accurately at the letters to play notes	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Pixel Paint

Choose an image, choose a colour, and get painting!

Skills: Exploration, Precision, Problem Solving



Creative

No Easy / Standard mode for this activity.

Skill	Objective	Not Observed	Developing	Achieved	Comments
Exploration	Looks at all areas of the picture	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Precision	Accurately targets remaining pixels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Problem solving	Picks correct colour for remaining pixels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	


Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Colour In

Choose a picture and colour it any way you like. Once you are done you can save a copy of your masterpiece.

Skills: Exploration, Choice making, Precision

In Easy mode: The images are less detailed, and the range of colours is reduced to give larger cells.



Creative

Skill	Objective	Not Observed	Developing	Achieved	Comments
Exploration	Looks at all areas of the picture	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Choice Making	Chooses colours to complete the picture	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Precision	Targets small parts of the picture	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

General	Low	Med	High	Comments	
Engagement (How interested is the user in the activity? Do they show enjoyment?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Independence (How much support do they need to access the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Success (How well do they succeed in the activity?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Name

Date

Facilitator

Time

Comments (e.g. Device position, user position, calibration)



Quick Record Sheet

Skills worked on:

Activity used

Observations

Activity used

Observations

Activity used

Observations

Activity used

Observations

Name

Date

Comments (e.g. Device position, user position, calibration)

Facilitator

Time

Comments (e.g. Device position, user position, calibration)

Find more Look Lab resources on
the Smartbox Hub:

hub.thinksmartbox.com/looklab



Smartbox

Document: Look Lab - Workbook - UK ENG - A
Language: English
Date of issue: June 2023

Copyright © 2023 Smartbox Assistive Technology

