

# Technical Specification

Remote Power Button 1 is an accessory included with Grid Pad 12 & 15. This allows you to turn your Grid Pad on or off.

Remote Power Button 2 is an accessory included with Grid Pad 10S. This allows you to turn your Grid Pad on or off. It also has switch access for switch scanning in Grid 3.



Name	Remote Power Button 1	Remote Power Button 2
Model	RPS-GP12	RPB2A
Price	Included with GP12, GP15	Included with GP10S
Warranty	Covered by Grid Pad warranty	Covered by Grid Pad warranty
Grid Pad Compatible	GP12, GP15	GP10S
Size	125 x 50 x 15 mm	125 x 50 x 15 mm
Weight	75g	75g
Batteries	Replaceable 2 x AAA	Replaceable 2 x AAA
Battery Life	2 years typical use	1 year typical use
Switch Port 1	Any press to turn device on Long hold to turn device off	Any press to turn device on Long hold to turn device off Switch Access in Grid 3  You can configure how this using Grid. Menu - Settings - Computer.
Switch Port 2	None	Switch Access in Grid 3
Other	Polycarbonate enclosure, with silicone button, drop tested to 1M	Polycarbonate enclosure, with silicone button, drop tested to 1M
Regulatory	CE RED FCC Part 15 b.	CE RED FCC Part 15 b.