

L 8


Slim

Specification

Front camera




Power / Sleep button



Volume controls

Mounting plate (optional /
sold separately)

Rear camera



Strap ring

Strap ring

Charging port

Contents

iPad Mini (8.3)
Smartbox Case
Strap
Bag
Screen protector
Charging cable
Mains adapter

To find out details about your iPad select

Settings - General - About

Battery life: 10 + hours
Access: Touch
Bluetooth switch (sold separately)
Weight: 500g
Mounting: Optional Rehadapt and Daessy (28g)
iPad: iPad Mini 6th generation

Document: L 8 - SL - Specification - US - A
Device: TS08E-L8
Language: English
Date of issue: September 2023

Copyright © 2023 Smartbox Assistive Technology