

Smartbox



Communication is life


Start here

Ready to go

Getting started with your device

This device is supplied **Ready to go**.


This means we have already prepared it to your requirements.



Turn on your device with the power button.

When you press the power button, the device will wake up.

There is no pass code on the device.



If your device has a speaker, power it on by pressing the connect button for 3 seconds.

Important things to know about your Device

Your device is pre-configured to meet your requirements and is managed by Smartbox Inc.

If your device has a speaker, it will be paired and ready to use.

We have installed any required apps, any addition information you may need is below.

Unlocking your device

Your device is managed by Smartbox Inc. If you would like to unlock your device to install other apps and services, you can contact us today.

Please have your device details to hand as we will use these to help you. These can be found in the information supplied with your device.

Online



hub.thinksmartbox.com/unlock

Get support

Phone: (844) 341-7386

Monday - Friday

Email: support@thinksmartbox.com